



SigNature[®]
molecular ink

V4221M



Unique to each customer, SigNature[®] molecular ink is created with the addition of a unique, secure molecular tag.



SigNature molecular ink is delivered via a secure SmartCartridge.[®]



SmartCartridge is paired with a Videojet 1860M continuous inkjet printer.



TAG products and packaging.



TEST to authenticate at any point along the supply chain.



TRACK to collect chain of custody data for traceability.



t +91 75060 01861 | marketing.india@videojet.com



t +1.631.240.8823 | bob.maddowell@adnas.com